

ADJUSTER[®] X

GLIDE
WEAR[®]


Available with GlideWear[®]
Technology. See pg. 25 for details.



*Shown with GlideWear[®]



8 ADJUSTABLE
CHAMBERS



*Standard Chambers Shown

E2622 / E2623**ADJUSTABLE SKIN PROTECTION****PRODUCT FEATURES**

- Designed with adjustable chambers, the Adjuster® X is the most versatile adjustable skin protection cushion on the market.
- Liner is made out of a breathable 3D mesh material and allows the full liner with air cells to be washed and dried.
- I.T. bony prominences have a dedicated chamber that can be partially or fully offloaded giving the maximum healing opportunity for any stage wound.
- Independently adjustable lateral chambers in the rear accommodate lateral pelvic asymmetries or loading the greater trochanters.
- Front chambers can be adjusted to handle seat surface height variances and alter the anti-thrust characteristics for better pelvic positioning.

CLINICIAN TIP

Versatile cushion that provides an ideal fit for a wide variety of body types

WEIGHT CAPACITY

20" Wide & Under	550 lbs
22" Wide & Wider	750 lbs

DIMENSIONS

Seat Thickness	4.5"
----------------	------

STANDARD SIZES

	10" W	12" W	14" W	16" W	18" W	20" W	22" W	24" W	26" W	28" W	30" W
10" D	10 x 10	12 x 10	14 x 10								
12" D	10 x 12	12 x 12	14 x 12	16 x 12							
14" D	10 x 14	12 x 14	14 x 14	16 x 14	18 x 14						
16" D	10 x 16	12 x 16	14 x 16	16 x 16	18 x 16	20 x 16	22 x 16	24 x 16			
18" D		12 x 18	14 x 18	16 x 18	18 x 18	20 x 18	22 x 18	24 x 18	26 x 18	28 x 18	30 x 18
20" D			14 x 20	16 x 20	18 x 20	20 x 20	22 x 20	24 x 20	26 x 20	28 x 20	30 x 20
22" D				16 x 22	18 x 22	20 x 22	22 x 22	24 x 22	26 x 22	28 x 22	30 x 22
24" D					18 x 24	20 x 24	22 x 24	24 x 24	26 x 24	28 x 24	30 x 24

Additional sizes available

ADJUSTER® X PART NUMBERS

Cover Style	Stretch-Air™	Comfort-Tek™
Standard	AJX-S-WWDD	AJX-F-WWDD
+ GlideWear®	AJX-SQ-WWDD	AJX-FQ-WWDD

WWDD to be replaced / **WW** = width & **DD** = depth in inches

ADJUSTER®**E2622 / E2623****ADJUSTABLE SKIN PROTECTION****PRODUCT FEATURES**

- Five fully adjustable chambers provide a custom fit with maintenance free Comfort Cells™.

CLINICIAN TIP

A high level of immersion provides a soft feel

DIMENSIONS

	Standard	Low
Seat Thickness	4"	3"

LOW SIZES

	16" W	18" W	20" W
16" D	16 x 16	18 x 16	20 x 16
18" D	16 x 18	18 x 18	20 x 18
20" D	16 x 20	18 x 20	20 x 20

STANDARD SIZES

	10" W	12" W	14" W	16" W	18" W	20" W	22" W	24" W
10" D	10 x 10	12 x 10	14 x 10					
12" D	10 x 12	12 x 12	14 x 12					
14" D	10 x 14	12 x 14	14 x 14	16 x 14				
16" D		12 x 16	14 x 16	16 x 16	18 x 16	20 x 16	22 x 16	24 x 16
18" D			14 x 18	16 x 18	18 x 18	20 x 18	22 x 18	24 x 18
20" D				16 x 20	18 x 20	20 x 20	22 x 20	24 x 20

Additional sizes available

Weight Capacity for All Sizes = 550 lbs

ADJUSTER® PART NUMBERS

Cushion Style	Thickness	Stretch-Air™	Comfort-Tek™
Standard	4"	AJ-S-WWDD	AJ-F-WWDD
Low	3"	AJ-SL-WWDD	AJ-FL-WWDD

WWDD to be replaced / **WW** = width & **DD** = depth in inches